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(EN)

The note below describes the steps to make a Radiologic Unit Replacement.

Manipulate covers with precaution to not stripe them (you can lay covers on soft area like towels for example).

WARNING The technician must comply with the standards currently in force while installing and using the radiology unit. Whenever necessary, he must use a shield or wear a lead-lined apron.

The device must be switched off and unplugged.

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
REMOVE THE COVERS / OUVERTURE DES CAPOTS

(EN)

Remove the adjustment knobs of the Frankfurt laser (a).
Remove the 4 screws and caps (b) while supporting the cover.
Remove the rear cover of the radiologic Unit.
Remove the 4 screws and caps (c) while supporting the rear cover.
Remove the rear cover.

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
REMOVE THE KODAK 8000 TUBEHEAD DEMONTAGE DE LA CUVE KODAK 8000

(EN)

If you don't have the cephalometric option, remove the primary collimator, attached to the radiologic unit by the 5 screws (b).
Disconnect the connector from the radiologic unit (c).
Remove the tubehead assembly (THA) bracket, 4 screws (e) and fit it on the new radiogenic unit.
Remove the defective radiogenic unit and replace it, securing it with the 4 screws (d).

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If you have the cephalometric option / Si vous avez un cephalostat

(EN)

Place firstly the Franckfurt laser horizontally.

Remove the pre-assembled collimator assembly with a 22mm flat spanner.


- (a) Remove the laser, 1 screw
- (b) Remove the primary collimator (to refit, note the right position of the window)
- (c) Remove the collimator detection board and its 2 spacers
- (d) Unplug the generator
- (e) Remove the THA bracket, 4 screws and fit it on the new radiogenic unit
- (f) Remove the defective radiogenic unit and replace it, securing it with the 4 screws

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TUBEHEAD REASSEMBLY / REMONTAGE DE LA CUVE

(EN)


Refit the components in the reverse order.

Check the dimension between the surface of the lower mounting bracket of the radiologic unit and the lower edge of the digital sensor. This dimension should be 461 mm. If this is not the case, insert washers as shims between the screws (d) in order to adjust the dimension.

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4

PROCEDURE OF REGENERATION OF THE TUBE AFTER A LONG PERIOD OF NON-USE


(EN)

If the product is new or not used for a long time, the X-ray tube may degas during this period (i.e. these gas molecules escape from materials constituting the tube). The phenomenon is particularly significant when the tube is not or little used. In the event of degasification, an immediate use with full power is likely to cause electrical arcing in the tube. To recreate the vacuum inside the tube, it is necessary to carry out the following procedure.

If the tube contains too much gas, there are risks of arcing that can be detected by electronics. So, following error messages may appear "E3", "E4", "E6", "E7", "E8" or "E10" In that case, wait for 15 minutes and start again the procedure from the beginning.

If after having carried out the procedure 3 times, an error message still appears, then the X-ray generator is defective.

During this procedure, it is recommended to measure the mA value. In case of fluctuation over 20% of the nominal value, start again the procedure from the beginning.

Reminder : mA calibration

1. Turn off the unit. Remove the top cover.
2. Remove J12 from the CJ664 (CPU board) and measure the resistance between pins 1 and 4 of the connector J12. The value must be between 244 and 253 Ohms.
3. Connect J12 back and connect a DC Voltmeter between the GND (PT15) and test point PT22 (RET – mA) on the CJ664 board.

THE VALUE MEASURED WILL BE OF 4mA PER VOLT DURING EXPOSURE

Last exposure carried out with tubehead less than 3 months ago :

Maintenance → **Calibration** → **Tubehead calibration**

- 60 kV – 2 mA – Make 3 exposures. Wait for 10 minutes
- 60 kV – 6.3 mA – Make 1 exposure. Wait for 1 minute
- 70 kV – 6.3 mA – Make 1 exposure. Wait for 1 minute
- 80 kV – 10 mA – Make 1 exposure. Wait for 1 minute
- 90 kV – 10 mA – Make 1 exposure

Last exposure carried out with tubehead more that 3 months ago :

Maintenance → **Calibration** → **Tubehead calibration**

- 60 kV – 2 mA – Make 10 exposures. Wait for 15 minutes
- 60 kV – 2 mA – Make 3 exposures. Wait for 1 minute
- 70 kV – 2 mA – Make 3 exposures. Wait for 1 minute
- 80 kV – 2 mA – Make 3 exposures. Wait for 1 minute
- 90 kV – 2 mA – Make 3 exposures. Wait for 15 minutes
- 60 kV – 15 mA – Make 3 exposures. Wait for 3 minutes between each exposure
- 70 kV – 15 mA – Make 3 exposures. Wait for 3 minutes between each exposure
- 80 kV – 15 mA – Make 3 exposures. Wait for 3 minutes between each exposure
- 90 kV – 15 mA – Make 3 exposures

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From the process interface, select :

**Maintenance
Calibration
mA Calibration**

Click on Start to Begin Calibration Process.

- A new window is displayed "mA Calibration- Step 1" (4mA position calibration). Make an exposure, note the displayed voltage, multiply this value by 4, and enter this value in the 'Milliammeter' field (this value must not vary by more than ± 1 mA from the set value. Otherwise, make another exposure or check the accuracy of your voltmeter). Click on Next.
- A new window is displayed: 'mA Calibration - Step 2' (12mA position calibration). Make an exposure, note the displayed voltage, multiply this value by 4, and enter this value in the 'Milliamperimeter' field (this value must not vary by more than ± 2 mA from the set value. Otherwise, make another exposure or check the accuracy of your voltmeter). Click on Next.
- The window indicates the new correction factors applied to the command and the mA measurement. The calibration process is complete. Click on Save.

If an error message is displayed during the calibration process, refer to the paragraph 'Error Messages' in your service guide.

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(EN)

From the process interface, select :

**Maintenance
Calibration
Heating Voltage Calibration**

- A window is displayed for calibrating the heater corresponding to each kV value. To do so, press the exposure button until the audible signal ceases. The microprocessor automatically measures the heater value for the first kV value, then moves on to the next value.
- Repeat the operation for the second value and so on for each kV value.
- When the last kV value has been calibrated, click on 'Save'.

If an error message is displayed during the calibration process, refer to the paragraph 'Error Messages' in your service guide.

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
(EN)

Place the panoramic unit in «service Mode»

1. Switch off the unit.
2. Press the 2  keys and, at the same time, switch on the device while holding down both keys for a few seconds.
The display should indicate 'Service Mode'.

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STEP 2

(EN)

Remove the cover from the digital sensor and place the scintillator in the secondary slot of the digital sensor (if necessary, use adhesive tape to hold it in position).

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(EN)

From the process interface, select :

Maintenance
Calibration
Tube Head Calibration

In the 'Parameters' window, select the radiological parameters required for correct illumination of the scintillator (75 kV, 10 mA for example).

IT IS RECOMMENDED THAT THE OPERATION SHOULD BE PERFORMED IN THE DARK IN ORDER TO CLEARLY DISPLAY THE FIELD ON THE SCINTILLATOR.

Make an exposure and check that the X-ray field is within the markings on the scintillator. If the adjustment is unsatisfactory, slightly loosen the 2 screws securing the primary collimator and adjust its position in order to rectify the fault. Make an exposure and check again.

Repeat the operation as often as necessary. When the adjustment is satisfactory, tighten the screws and refit the covers.

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
Aligning the beam

(EN)

Place a metal rod instead of the biteblock rod on the pan chinrest.

Move the arm to the median position, so that the center line of the digital sensor, the metal rod and the center line of the radiogenic unit are in alignment.

Immobilize the arm in this position, using the screw (a) located beneath the upper cover. This screw passes through the CJ653 board.

To make adjustments using the scintillator, select 'Check Generator' mode.

To make adjustments without using the scintillator, select 'Check X-ray Beam' mode. See above.

Make an exposure and check that the shadow of the metal rod is centered on the scintillator or on the image.

If the adjustment is unsatisfactory, remove the lower cover from the control panel.

Slightly loosen the 2 screws (a) securing the chinrest assembly and adjust its position in order to rectify the fault. Make an exposure and check that the metal rod is centered.

Repeat the operation as often as necessary. When the adjustment is satisfactory, tighten the screws and refit the cover.

To exit from the 'Check Generator' or 'Check X-ray Beam' modes, select a type of image in the Program window.

Do not forget to release the arm by removing the screw (a) passing through the CJ653 board.

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
(EN)

Remove the bite block and the rod.

Position the ball phantom film and check that the graduated scale on the chinrest is set to the '0' position.

Make a panoramic exposure at 60 kV 2 mA. The image should be symmetrical and the balls should be circular.

If the balls are not circular, move the graduated scale at millimeter intervals, in one direction or another, making an exposure each time, until the central ball is circular.

Mark the position of the scale at which the center ball is circular and remove the chinrest.

Move the label ref. JP720 (affixed to the chinrest assembly) so that it meets this position.

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
Adjusting the lasers

(EN)

- Switch on the lasers. Check that the sagittal median and canine lasers intersect perpendicularly over the center ball.
- Slightly loosen the screws retaining the laser mounting (a) *canine* or (b) *sagittal median*, then adjust the laser so that the beam is correctly aligned on the center ball. Tighten the screw, when it is correct.
- Check that the laser of the Frankfurt plane is perpendicular to the center line of the digital sensor.
- Slightly loosen the screw retaining the laser (c) to be adjusted, then adjust the laser so that the beam is perpendicular to the vertical center line of the digital sensor. When it is correct, tighten the screw.

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